NECK TRACTION

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**Unbelievable Superior Results of Pain Relieve
within a Minutes By Using an AIR NECK TRACTION...**

Neck traction devices may improve the blood supply to neck muscles, tendons, and ligaments and promote relaxation by relieving pressure, stiffness and stress in the neck, head, and shoulders.

Tight muscles can restrict blood flow causing a release of inflammatory chemicals which create further irritation, tension and pain.

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NECK TRACTION

Forget your old ways of Neck Traction ... and use a very simple ways to get the same ultimate result of neck pain relive and help to relax your muscular Spasms with increases joint mobility of the cervical spine at any time and at Any Place (Home, Office, Car, Etc...)

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
New Ways with Air Neck Traction ...

Pain Relieve, Pinched Neck Nerve ...

Qualities are Guaranteed
www.remalcare.com
Made In China

Why Neck Traction

Neck traction may be helpful for capsulitis, cervical facet joint syndrome, and in helping to relax muscular spasms.

Traction may also benefit patients with osteoarthritis by promoting increased movement of the cervical spine and decreased pressure on the facet joints.

Cervical traction increases joint mobility of the cervical spine for patients with hypomobility. Degenerative disc disease often leads to radiculopathy and cervical pain.

Air Neck Traction alleviates pressure on the diseased discs to promote decreased pressure on peripheral nerve roots and enhance nutrient delivery to the disc.

Instructions For Using The Air Neck Traction

You will achieve the best results with regular applications. Begin your neck traction program with short traction sessions (3-5 minutes each).

After the first week, you can gradually increase the amount of pressure and length of the traction sessions by a few minutes. Sessions should always be comfortable.

If you notice any increase soreness and symptoms returning, decrease the length of your session.

Take things slowly in the beginning for optimum results.