

DM-50 - Home Use / Hospital / Clinic UV Phototherapy Lamps-Psoriasis,Vitiligo

Home Use UV Phototherapy Lamps-Psoriasis,Vitiligo. Practically New! Only had about 2 hours of use (max). No longer need it because my condition has gone into remission and no longer need it.

The Lamp was supplied as standard with UVB-Narrow Band 311nm bulbs but was also available with either UVA or UVB from suppliers website.

The Lamp is easy to use and also ideal for in-home phototherapy. The standard Narrow Band UVB-311(Philips) is best suited for psoriasis, vitiligo, eczema, acne, rosacea, alopecia hair loss and similar conditions.



1. No prescription or doctor's letter required
2. The Lamp is a light weight, ergonomic, hand held UV phototherapy lamp. Perfect for treatment of Psoriasis, Vitiligo, Dermatitis, Acne or Alopecia hair loss.

Features:

1. no prescription or doctor's letter required to order.
2. Phillips UVB NB-311 narrow band bulbs are standard.
3. comb attachment included for scalp treatment.
4. uniyt instructions, timer and UV safety goggles include.
5. bulb life at least 300 hours or several years.

Technical Data:

Mains voltage	: 210-230V /50Hz~60Hz
Power consumption	: 0.17A/cosφ0.31
Connecting cable	: 3 m with plug ballast unit
Built-in fuses	: 2 micro fuses T6.3H 250V
Radiant tube	: 1×PL-S 9W/01 lamp
Protective tube	: highly impact-resistant plastics, transparent
Net weight	: 0.68kg(without packaging)
Packed weight	: 2kg (with packaging)
Packed dimensions	: 36(L) x 30(W) x 10(H)cm



**The Narrow Band - Lamp - UVB-311(Philips) Made In Poland
Handel Made In China**

**DM-50 - Home Use / Hospital / Clinic
UV Phototherapy Lamps-Psoriasis,Vitiligo**

Special characteristics:

Light-guiding, rounded comb attachment for a smooth and precise spot treatment of the scalp. Wide treatment area. Lightweight and handy construction. Short treatment time. Easy and quick to clean. Ideal for in-home therapy. High-quality workmanship. Ageless design Therapy!

What Skin Diseases May Be Treated By Therapy Comb?

You should always agree the use of the therapy comb with your dermatologist. He decides if in your case a treatment is useful or not.

Mainly the therapy comb is used for treatment of Psoriasis of the scalp. But also in Seborrheic Eczema and other stubborn dandruffs (Pityriasis simplex) the therapy comb may be used successful.

After removing the comb the device may be used as therapy pencil.

This is mainly a possibility for treatment of skin areas that are not sufficient irradiated during normal whole body phototherapy.

Such skin areas are e. G. Armpit, genital area or anus.

You have to note that the device must not be used without moving because otherwise skin irritations may occur. It is practical to move the device in this kind of treatment too.

The therapy device is partly useful for treatment of psoriasis around the fingernails. Mainly the environment of the nail can be positively influenced but not the fingernail itself.

For treatment of the fingernail diseases the therapy pencil may be used for local photo chemotherapy (PUVA) too. This treatment must always be done in a dermatologist's practise. For use in PUVA treatment a special manual exists but it is only passed on to dermatologists.

About any other application only a dermatologist may decide.

Note:

1. In narrowband-UVB-therapy in every treatment the treatment should be extended as long as just still no sunburn (erythema) is caused.
2. Before treating Psoriasis remove as much flaking skin as possible because the radiation hardly penetrates the flaking skin.
3. For treatment of not hairy skin without comb the device should be moved in a minimum distance to the skin (roughly 1 centimetre).



Eczema Treatment Phototherapy

علاج الأكزيما الجلدية بالضوء



Skin Disease & Atopic Dermatitis Treatment

علاج أمراض و التهاب الجلد الاستشرائية



Ultraviolet Lamp for Psoriasis

الأشعة فوق البنفسجية لعلاج الصدفية



Vitiligo Phototherapy

علاج البهاق بالضوء

